

Klenow Reaction and Cy Dye coupling for *P. falciparum* genomic DNA
9/29/2003

3 µg genomicDNA, sheared to between 500 bp-1.5 kb fragments (2 µg genomicDNA/hybe)
 50 µg pdN9 primer
 1x EcoPol buffer (NEB)
 1x dNTP mix (75 µM dATP, 15 µM dTTP, 60 µM aa-dUTP, 15 µM dGTP, 15 µM dCTP)
 50 U klenow/50 µl reaction

	Volume (1 rxn)	Volume (5 rxns)	Volume (9 rxns)	Volume (11 rxns)
200 ng/µl gDNA	15 µl			
10 mg/ml pdN9	0.25 µl			
Total volume	15.25 µl			
10x EcoPol buff	5 µl	25 µl	45 µl	55 µl
10x dNTP mix	5 µl	25 µl	45 µl	55 µl
Klenow	1 µl	5 µl	9 µl	11 µl
ddH ₂ O	23.75 µl	118.75 µl	213.75 µl	316.25 µl
Total volume	34.75	173.75	312.75	437.25
Final vol/rxn	50 µl	50 µl	50 µl	50 µl

1. Boil gDNA/pdN9 for 5 minutes at 99°C, then cool on ice for 5 minutes.
2. Add 2nd mixture. Mix well.
3. Incubate: 37°C, 2-3 hrs.
4. Add additional 1 µl klenow, leave at 37°C, 2-3 hrs.

20x dNTPs (1 dGTP:1 dCTP:5 dTTP/aa-dUTP:5 dATP):

<i>100 mM stock solutions</i>	<i>20x final conc.</i>
dATP 3.75 µl	1.5 mM
dTTP 0.75 µl	0.3 mM
aa-dUTP 3.00 µl	1.2 mM
dGTP 0.75 µl	0.3 mM
dCTP 0.75 µl	0.3 mM
<u>ddH₂O 241 µl</u>	
total vol. 250 µl	

10x dNTPs (1 dGTP:1 dCTP:5 dTTP/aa-dUTP:5 dATP):

<i>100 mM stock solutions</i>	<i>10x final conc.</i>
dATP 1.875 µl	0.75 mM
dTTP 0.375 µl	0.15 mM
aa-dUTP 1.50 µl	0.60 mM
dGTP 0.375 µl	0.15 mM
dCTP 0.375 µl	0.15 mM
<u>ddH₂O 245.5 µl</u>	
total vol. 250 µl	

Clean up Klenow reactions with Zymo columns (200 µl DNA binding buffer/50 µl rxn mix). Wash 2x with 200 µl wash buffer and elute with 12 µl ddH₂O. Add 1 µl 1M NaBicarb, pH 9.0 (final conc. 0.1 M NaBicarb). Couple DNA to appropriate Cy dye 2 hrs to O/N in the dark at RT. Clean up coupled Klenow reactions with Zymo columns (200 µl DNA binding buffer/10 µl mix as before, elute with 16 µl EB buffer (Qiagen)). Ready for hybridization to DNA microarrays.

40 µl hybe mix (for thick lifter slips):

Boil mixture 2 mins, 99°C.
 Cool by spinning at RT.
 Apply entire mix to array
 Hybridize at 65°C, O/N

30 µl labeled DNA (15 of each-cy3 & cy5)

6 µl 20x SSC
 2 µl 10% SDS
1 µl 1M Hepes, pH 7
 40 µl total volume